

## Two-Tone Disclosing Tablets



Without  
Erythrosine  
and  
Gluten




- For quick and easy identification
- Dual-colour detection: tablets turn **newer plaque red** and **older plaque blue** and help showing the areas consistently missed while brushing.
- Ideal for use at home. Help educate children on how to improve their at-home oral care routines.
- Fruity taste

### Instructions for use:

1. Chew one tablet.
2. Wipe teeth with tongue.
3. Rinse mouth. Spit out.
4. Brush and floss teeth thoroughly to remove all stain.

<b>Item #</b>	<b>Disclosing Tablets</b>
295718	individually wrapped, 250 ct. box



## Paediatric Bibs



Five fun-themed characters



### Get into character with Zooby bibs!

Each Zooby® bib features the body of an animal and acts as costume that your paediatric patient will love!

- Characters: Penelope - the Penguin, Tobi - the Turtle, Francesca - the Flamingo, Talon - the Tiger and Allie - the Alligator
- Paediatric size (24 x 23 cm)
- 1-ply tissue and 1-ply poly



**Item #** Paediatric Bibs - 100 ct.  
295200 Assorted pack, 5 designs

ZB-GB-TB/O921