

# YOUNG™

## DISPOSABLE PROPHY ANGLES

Young prophy angles are designed to access hard-to-reach areas with ease, while a smooth, round head design make polishing enjoyable for both hygienists and patients.

Reduced risk of infection, no damage to the costly handpiece due to abrasive pastes. They save time and money and are single use - no disinfection required.

17° contra angle offers improved ergonomics



Classic Petite Web™  
295683

Contra Petite Web™  
295682

Contra Elite™  
295680

Classic Elite™  
295681

**Slim Neck**  
Diameter for  
Greater Visibility

Enhanced  
**Grip Features**

**Rounded Head**  
for Improved  
Maneuverability

Sleek Design for  
**Clinician and  
Patient Comfort**

All Young prophy cups are made with a proprietary rubber. FREE FROM LATEX AND BPA.

### Comparison old and new generation Disposable Propy Angles



25% Slimmer Neck Design



Sleek Design with Reduced  
Surface Area



58% Greater Clearance with  
rounded Head Design

#### Item # Disposable Propy Angle - 125ct

- 295680 Contra Elite™, Soft, Latex Free
- 295681 Classic Elite™, Soft, Latex Free
- 295682 Contra Petite Web™, Soft, Latex Free
- 295683 Classic Petite Web™, Soft, Latex Free

# YOUNG PROXEO HANDPIECE

The Young Proxexo handpiece for Young disposable prophylaxis angles allows for user-friendly and low maintenance application.

LIGHTWEIGHT  
ADAPTER  
MINIMIZES  
FATIGUE

User-friendly **automatic clamping**  
chuck system

Doriot Shank System

4:1 conversion for any  
E-Style electric motor

For use with  
DPA 's with  
plastic shafts

**Sterilisable** and  
for use in **thermo**  
washer disinfectant



Young Proxexo Handpiece HP-44M



HP-44M with Classic Petite Web™



HP-44M with Contra Elite™

296063/0219