

# YOUNG®

## DISPOSABLE PROPHY ANGLES

Young prophy angles are designed to access hard-to-reach areas with ease, while a smooth, round head design make polishing enjoyable for both clinicians and patients.

Reduced risk of infection, no damage to the handpiece due to abrasive pastes.

They save time and money and are single use - no disinfection required.

Proprietary rubber  
free from latex  
and BPA

17° contra angle for  
improved ergonomics



Rounded Head  
for improved  
manoeuvrability

Slim Neck  
diameter offers  
greater visibility

Sleek design for  
clinician and  
patient comfort

Enhanced  
grip features

295683  
Classic Petite  
Web™

295682  
Contra Petite  
Web™

295680  
Contra Elite™

295681  
Classic Elite™

### Comparison old and new generation disposable prophy angles



25% slimmer neck design



Sleek design with reduced  
surface area



58% greater clearance with  
rounded head design

### Item # Disposable Prophy Angles - 125 ct.

- 295680 Contra Elite™, soft, without Latex and BPA
- 295681 Classic Elite™, soft, without Latex and BPA
- 295682 Contra Petite Web™, soft, without Latex and BPA
- 295683 Classic Petite Web™, soft, without Latex and BPA

# YOUNG PROXEO HANDPIECE

The Young Proxexo handpiece for Young disposable prophy angles allows for user-friendly and low maintenance application.



Lightweight adapter minimises fatigue

User-friendly **automatic clamping chuck system**

**Doriot** shank system

4:1 conversion for any E-Style electric motor

For use with **DPA 's with plastic shafts**

**Sterilisable** and for use in **thermo washer disinfectant**

## ALSO AVAILABLE AS STARTER KIT

375 Disposable Prophy Angles  
+ 1 Young Proxexo Handpiece



Item #	Description Starter Kit
YDNT01	Starter Kit YDNT01, Contra Elite™
YDNT02	Starter Kit YDNT02, Classic Elite™
YDNT03	Starter Kit YDNT03, Contra Petite Web™
YDNT04	Starter Kit YDNT04, Classic Petite Web™

**Risk-free trial!**  
Disposable Prophy Angles  
and Proxexo handpiece.



YD-GB-DPA/0921