

YOUNG

INNOVATIONS

PRODUCT CATALOGUE

Europe

ORAL HYGIENE

PEDIATRIC

RESTORATIVE

ABOUT YOUNG INNOVATIONS

Young Innovations continually seeks to equip oral health professionals with the knowledge, products and services they require to meet the needs of patients — ultimately our shared and most important customers.



Since the early 1900s, Young Innovations has been delivering premium products that continue to push the dental industry forward. We have spent the last century cultivating a strong connection with the dental community and developing a deep understanding of what clinicians need as the industry has evolved. Dental professionals have always been our primary focus and, as such, their input directly influences our product offerings in the preventive, restorative, endodontic and diagnostic segments. Our consistent focus on our customers and on product innovation has led us to become a leading designer and manufacturer of preventive dental products, restorative accessories and instruments dental professionals can rely on to enhance their clinical success.

YOUNG
INNOVATIONS



TABLE OF CONTENTS



ORAL HYGIENE 4-9

Prophy Cups4-5

Young™ Prophy Cups

Disposable Prophy Angles6-7

Young™ Disposable Prophy Angles

Young™ Original Disposable Prophy Angles

Miscellaneous 8-9

2 Tone™ / Trace® Disclosing Agents

Surg-O-Vac™ / Vent-O-Vac™ Moisture Control



PEDIATRIC 10-11

Disposable Prophy Angles 11

Zooby® Disposable Prophy Angles

Prophy Disclothing Tablets 11

Zooby® Disclothing Tablets

Prophy Bibs 11

Zooby® Bibs



RESTORATIVE 12-17

Applicators 12-15

Microbrush® Plus Dispenser Kit

Microbrush® Plus Dispenser Refill

Microbrush® Tube Series

Microbrush® X Applicator

Micro-Stix™ Adhesive Tip Applicator

TRU™ by Microbrush® Applicator

UltraBrush® Bristle Brush

Moisture Control 16-17

DryTips®

Reflective DryTips®

Dri-Aid™

Dri-Aid™ Silver



INSTRUMENT CASSETTES 18

GALAXIE Instrument cassettes 18

YOUNG™

ELITE™ PROPHY CUPS

Young offers dental professionals the widest variety of prophy cups in both latch and screw type to suit every need. Young guarantees that each cup will provide the same durable, consistent and reliable performance for each polishing procedure. Cups are color coded for simplicity.


All Prophylaxis
Cups are Made
With Proprietary
**LATEX AND BPA-
FREE**
Rubber!

Proprietary rubber offers
greater flare for a better
polishing experience



Internal webbing conforms
to the tooth surface and
reduces splatter

Exterior ridges remove
interproximal stains



Engineered specifically for the
unique requirements of the
European dental professional



To learn more about our products, visit www.youngdental.eu

PROPHY CUPS

Young™ Prophylaxis Cups

- » Each cup is designed to be used on any patient but has specific features that make it more useful in different clinical situations (mouth sizes, stain level, periodontal health status, etc.)
- » We have unique dual-colored cups for easy identification. Our outside cup color represents firmness level (● Blue=Firm, ● Purple=Soft, ● Pink=Extra Soft)
- » All prophylaxis cups are made with a proprietary latex-free rubber
- » Our prophylaxis cups also come individually sealed in unique blister packaging. This greatly reduces the risk of cross contamination, and reinforces the single-use concept

Cup Type Overview

Product	Name	Item#	Quantity	Attachment	Latex-Free	Firmness	Size	Exterior Ridges
	Young™ Web	295661	40ct	Screw	✓	Firm	Long	
	Young™ Web	295662	40ct	Screw	✓	Soft	Long	
	Young™ Turbo Plus™	295663	40ct	Screw	✓	Firm	Long	
	Young™ Turbo Plus™	295664	40ct	Screw	✓	Soft	Long	
	Young™ Elite™	295665	40ct	Screw	✓	Soft	Short	✓
	Young™ Elite Flex™	295666	40ct	Screw	✓	Extra Soft	Short	✓
	Young™ Petite Web™	295667	40ct	Screw	✓	Soft	Short	
	Young™ Petite Flex	295668	40ct	Screw	✓	Extra Soft	Short	
	Young™ Pointed Polisher	295669	40ct	Screw	✓	Soft	Long	✓
	Young™ Web, Long	295671	40ct	Latch/RA	✓	Firm	Long	
	Young™ Web, Long	295672	40ct	Latch/RA	✓	Soft	Long	
	Young™ Turbo Plus™	295673	40ct	Latch/RA	✓	Firm	Long	
	Young™ Turbo Plus™	295674	40ct	Latch/RA	✓	Soft	Long	
	Young™ Elite™	295675	40ct	Latch/RA	✓	Soft	Short	✓
	Young™ Elite Flex™	295676	40ct	Latch/RA	✓	Extra Soft	Short	✓
	Young™ Petite Web™	295677	40ct	Latch/RA	✓	Soft	Short	
	Young™ Petite Flex	295678	40ct	Latch/RA	✓	Extra Soft	Short	

YOUNG™

DISPOSABLE PROPHY ANGLES

Young prophy angles are designed to access hard-to-reach areas with ease, while a smooth, round head design make polishing enjoyable for both hygienists and patients. Disposable prophy angles reduce the risk of cross-contamination and save time because they are single use and don't require disinfection.



Save Time and Money with a DPA by reducing processing cost and wear of handpiece.



Rounded Head
for Improved
Manoeuvrability

Sleek Design for
Clinician and
Patient Comfort

Enhanced
Grip Features



Petite Web™

Contra Elite™

Slim Neck
Diameter for
Greater Visibility

**17 degree
Contra
Angle**
Offers
Improved
Ergonomics



Lightweight adapter minimizes fatigue and accommodates Young disposable prophy angles
4:1 conversion for any E-style electric motor

To learn more about our products, visit www.youngdental.eu

DISPOSABLE PROPHY ANGLES

Young™ Disposable Propphy Angles

- » 25% slimmer neck design allows for easy reach and access
- » Designed to reach different areas of the mouth with ease, enhancing clinician and patient comfort
- » The Contra body style has an ergonomically designed 17° angulation to reduce muscle fatigue and improve hand positioning
- » Dual-colored cups easily identify Young prophylaxis angles



Description	Quantity	Item#
Contra Elite™, Soft, Latex Free	125ct	295680
Classic Elite™, Soft, Latex Free	125ct	295681
Contra Petite Web™, Soft, Latex Free	125ct	295682
Classic Petite Web™, Soft, Latex Free	125ct	295683

Young™ Original Disposable Propphy Angles

- » We have developed an entire portfolio of disposable angles, providing dental professionals with options to address their every need
- » Our angles are specially constructed with an embedded screw to deliver a smooth, concentric rotation. This proprietary design keeps the cup in close contact with the tooth for a dependable polish every time
- » Each angle is lubricated for a reliable, smooth, vibration-free operation
- » These angles are uniquely designed to fit comfortably and securely in the hand and feature Young's high quality prophylaxis cups



Description	Quantity	Item#
Young™ Contra Turbo Plus™ LF - Soft	200ct	153620
Young™ Classic Brush - Firm	100ct	130110
Young™ Contra Petite Web™ LF - Soft	200ct	154620
Young™ Classic Traditional Web™ - Soft	200ct	131120
Young™ Contra Elite Original™ LF - Soft	200ct	155620
Young™ Classic Traditional Web™ - Firm	200ct	131320
Young™ Classic Elite Original™ LF - Soft	1200ct	135612
Young™ Contra Turbo Plus™ - Soft	200ct	153120
Young™ Classic Turbo™ - Soft	200ct	132320
Young™ Classic Elite Original™ LF - Soft	200ct	135620
Young™ Contra Turbo Plus™ LF - Firm	200ct	153720
Young™ Classic Junior® - Soft	200ct	131420
Young™ Contra Turbo Plus™ - Firm	200ct	153320
Young™ Classic Traditional Web™ LF - Soft	200ct	131620
Young™ Classic Elite Flex™ - Soft	200ct	138920
Young™ Classic Petite Web™ LF - Soft	200ct	134620

Young™ Disposable Pointed Polisher and Brush

- » Pointed Polisher has a unique pointed tip to remove tough stains
- » Pointed Polisher is flexible and soft enough to polish around brackets and under orthodontic wires
- » Brush removes stain from occlusal pits and fissures
- » Brush is great for preparing tooth surfaces prior to sealant and orthodontic procedures



Description	Quantity	Item#
Young™ Classic Pointed Polisher LF	100ct	137710
Young™ Classic Brush	100ct	130110



PEDIATRIC DISPOSABLE PROPHY ANGLES

Introduce children to Penelope the Penguin, Francesca the Flamingo, Chase the Cheetah, Talon the Tiger and Zoe the Zebra with our fun, animal-themed prophy angles! These prophy angle characters can be used to engage young patients, helping transform an ordinary dental visit into a more captivating encounter where pediatric patients don't just open wide, they open wide with a roar!



Putting the fun in pediatric dental care!



To learn more about our products, visit www.youngdental.eu

PEDIATRIC

Zooby® Disposable Propphy Angles

- » Introduce children to their favorite characters with our fun, animal-themed prophylaxis angles
- » These personified prophylaxis angles can be used to engage young patients, helping transform an ordinary dental visit into a captivating encounter where pediatric patients don't just open wide, they open wide with a roar!



Box includes
assortment of
all 5 animals!

Description	Quantity	Item#
Assorted pack, 5 animal designs	100ct	575010
Assorted pack, 5 animal designs	500ct	575050

Zooby® Pediatric Bibs

- » Get into character with Zooby Bibs!
- » Each Zooby pediatric bib features the body of an animal and acts as a costume that your pediatric patient will love!
- » Non-Permeable
- » Pediatric size (24 x 23 cm)
- » Five fun-themed Zooby characters: Penelope - the Penguin, Zoe - the Zebra, Francesca - the Flamingo, Chase - the Cheetah and Talon - the Tiger



Description	Quantity	Item#
Assorted pack, 5 animal designs, 20 of each	100ct	295200

Zooby® Disclosing Tablets

- » Advanced multicolored disclosing agents that allow quick and easy identification
- » Zooby® Disclosing Tablets turn newer plaque red and older plaque blue, and help educating children on how to improve their at-home oral care routines
- » Gluten-free and without erythrosine (E127). Fruity taste.



WITHOUT
ERYTHROSINE
AND GLUTEN

Description	Quantity	Item#
Plaque control, dual-color detection, individually wrapped	250ct	295718



Microbrush[®] International



APPLICATORS

The award-winning Microbrush[®] applicator is the top choice for precise application in areas with limited access. The stiff head allows for scrubbing and the strong bending portion provides complete control.

Available in multiple colors to identify multi-step procedures, these applicators can be used with cements, etchants, sealants, bonding agents, hemostatic solutions and more.



REGULAR — 2.0 mm



FINE — 1.5 mm



SUPERFINE
1.0 mm



ULTRAFINE
0.5 mm



APPLICATORS

Microbrush® Plus Dispenser Kit

- » Stiff head to allow for scrubbing
- » Strong bending portion for improved control
- » Multiple applicator colors available to identify different steps in procedure
- » Sion™ dispenser with built-in antimicrobial properties, specially designed to release one applicator at a time
- » Offered in 4 Sizes: Regular, Fine, Superfine and ultrafine



Description	Quantity	Item#	Color
Microbrush® Plus Dispenser Kit Fine	50ct	MPD	Yellow
Microbrush® Plus Dispenser Kit Regular	400ct	PR400KIT	Assorted
Microbrush® Plus Dispenser Kit Fine	400ct	PF400KIT	Assorted
Microbrush® Plus Dispenser Kit Superfine	400ct	PSF400KIT	White

Microbrush® Plus Dispenser Refill

- » Stiff head to allow for scrubbing
- » Strong bending portion for improved control
- » Multiple applicator colors available to identify different steps in procedure
- » 100 applicators per tray; 4 trays per pack
- » Offered in 4 Sizes: Regular, Fine, Superfine, Ultrafine



Description	Quantity	Item#	Color
Microbrush® Plus Dispenser Refill Regular	400ct	PR400	Assorted
Microbrush® Plus Dispenser Refill Regular	400ct	PR400BL	Blue
Microbrush® Plus Dispenser Refill Regular	400ct	PR400GR	Green
Microbrush® Plus Dispenser Refill Regular	400ct	PR400PE	Peach
Microbrush® Plus Dispenser Refill Regular	400ct	PR400PU	Purple
Microbrush® Plus Dispenser Refill Fine	400ct	PF400	Assorted
Microbrush® Plus Dispenser Refill Fine	400ct	PF400PI	Pink
Microbrush® Plus Dispenser Refill Fine	400ct	PF400YE	Yellow
Microbrush® Plus Dispenser Refill Superfine	400ct	PSF400	White
Microbrush® Plus Dispenser Refill Ultrafine	400c	PU400TE	Teal
Microbrush® Plus Dispenser Refill Ultrafine	400c	PU400OR	Orange

APPLICATORS

Microbrush® Tube Series

- » Stiff head to allow for scrubbing
- » Strong bending portion for improved control
- » Multiple applicator colors available to identify different steps in procedure
- » 100 applicators per dispensing tube; 4 tubes per pack
- » Offered in 4 Sizes: Regular, Fine, Superfine, Ultrafine



Description	Quantity	Item#	Color
Microbrush® Tube Regular	400ct	MRA400	Assorted
Microbrush® Tube Regular	400ct	MRB400	Blue
Microbrush® Tube Regular	400ct	MRG400	Green
Microbrush® Tube Regular	400ct	MRP400	Purple
Microbrush® Tube Regular	400ct	MRS400	Peach
Microbrush® Tube Fine	400ct	MFA400	Assorted
Microbrush® Tube Fine	400ct	MFP400	Pink
Microbrush® Tube Fine	400ct	MFY400	Yellow
Microbrush® Tube Superfine	400ct	MSF400	White
Microbrush® Tube Ultrafine	400ct	MUT400	Teal
Microbrush® Tube Ultrafine	400ct	MUO400	Orange

Microbrush® X Applicator

- » Extended and extra thin for use in confined spaces
- » Ideal for endodontic treatment
- » Designed with stiffer fibers for scrubbing
- » Designed for even placement of material to inner surface of tooth wall
- » Bending portion bendable up to 90 degrees
- » Refill trays fit the MB Sion dispenser – Kit contains 1 sion dispenser



Description	Quantity	Item#	Color
Microbrush® X Dispenser Kit	100ct	MPX	Black
Microbrush® X Refill	100ct	PX100	Black

Micro-Stix™ Adhesive Tip Applicator

- » Adhesive tip applicators designed to handle delicate items
- » Strong adhesive that securely holds larger items
- » Designed to leave no residue on the tooth



Description	Quantity	Item#	Color
Micro-Stix™	64ct	STIX64B	Blue

APPLICATORS

TRU™ by Microbrush® Applicator

- » 12-centimeter handle with stay-put bending portion for control
- » Advanced adhesive technology designed for greater consistency
- » Double-bending tip
- » New Ultrafine 0.5-millimeter head size designed for greater access
- » Sold in individual dispensing packs of 80 ct applicators



Description	Quantity	Item#	Color
TRU™ Microbrush® Tube Regular (2.0mm)	80ct	291159	Blue
TRU™ Microbrush® Tube Regular (2.0mm)	80ct	291160	Green
TRU™ Microbrush® Tube Fine (1.5mm)	80ct	291161	Yellow
TRU™ Microbrush® Tube Fine (1.5mm)	80ct	291162	Pink
TRU™ Microbrush® Tube Superfine (1.0mm)	80ct	291163	Gold
TRU™ Microbrush® Tube Superfine (1.0mm)	80ct	291164	Silver
TRU™ Microbrush® Tube Ultrafine (0.5mm)	80ct	291166	Orange

UltraBrush® Bristle Brush

- » Bristle brush applicators for spreading and smoothing materials
- » Bends close to tip for improved control
- » Chemical-resistant fibers designed to hold up to the harshest dental materials
- » Available in two sizes: Regular size for better spreading action and the fine size for more precise applications
- » Refill trays fit the Microbrush Sion dispenser
- » Kit contains 1 sion dispenser



Description	Quantity	Item#	Color
UltraBrush® Dispenser Kit Regular	100ct	U2D	Green
UltraBrush® Refill Regular	200ct	U2R200	Gold / Green
UltraBrush® Dispenser Kit Fine	100ct	U1D	Blue
UltraBrush® Refill Fine	200ct	U1R200	Blue / Orange



Microbrush[®] International

MOISTURE CONTROL

DryTips[®] provide ultimate saliva absorbency and don't lose moisture when handled. They far outlast cotton rolls and comfortably stay in place to provide complete coverage of the parotid gland and buccal mucosa.



To learn more about our products, visit www.microbrush.eu

MOISTURE CONTROL

DryTips®

- » Effective coverage of parotid orifice and buccal mucosa
- » Polyethylene film for cutting protection
- » Enables easier access to molar region
- » Stays comfortably in place, yet is easily removed with water spray
- » Flexible, comfortably accommodates to cheek movements
- » Ideal for impressions, occlusal adjustments, orthodontics, composites, fissure sealing and more



See instructional video at www.microbrush.eu

Description	Quantity	Item#	Color
DryTips® Small	50ct	291543	White
DryTips® Large	50ct	291542	White

Reflective DryTips®

- » Effective coverage of parotid orifice and buccal mucosa
- » Reflective coating for increased visibility and to aid in retraction
- » Polyethylene film for cutting protection
- » Enables easier access to molar region
- » Stays comfortably in place, yet is easily removed with water spray
- » Flexible, comfortably accommodates to cheek movements
- » Ideal for impressions, occlusal adjustments, orthodontics, composites, fissure sealing and more

NEW!



See instructional video at www.microbrush.eu

Description	Quantity	Item#	Color
Reflective DryTips® Small	50ct	291727	Silver
Reflective DryTips® Large	50ct	291726	Silver

Dri-Aid™

- » Thin, absorbent wafers
- » Allows for better access and keeps the working area dry



Description	Quantity	Item#	Color
Dri-Aid™ Small	250ct	331125	White
Dri-Aid™ Small	750ct	331190	White
Dri-Aid™ Large	250ct	331320	White
Dri-Aid™ Large	750ct	331375	White

Silver Dri-Aid™

- » Silver Dri-Aid™ coating reflects light for improved visibility
- » Thin, absorbent wafers
- » Allows for better access and keeps the working area dry



Description	Quantity	Item#	Color
Silver Dri-Aid™ Small	250ct	332125	Silver
Silver Dri-Aid™ Small	750ct	332190	Silver
Silver Dri-Aid™ Large	250ct	332325	Silver
Silver Dri-Aid™ Large	750ct	332375	Silver

INSTRUMENT CASSETTES

GALAXIE Cassettes

The gold-standard solution for instrument reprocessing.

A new range of cassettes - ergonomic and smart - offers an optimal solution for the classification, cleaning and sterilization of dental hand-instruments, while increasing safety and efficiency in the practice.



DIN cassette for 20 instruments
PP40373-3

High Quality Materials:

- » Electropolished stainless steel makes the cassette edges smooth and protects against corrosion.
- » Only food grade silicone used.

Innovative and ergonomic closing latch:

- » Safe and easy-to-use locking system, allows one-handed opening.

New design for a better cleaning process:

- » The new and innovative laser cutting pattern guarantees maximum efficiency during the cleaning and sterilization process whilst ensuring stability and strength of the cassette.
- » Open lattice design specially designed for washer disinfectant use, enabling maximum water through flow.
- » Smooth round corners and slotted edges increase drainage and reduce drying time.

Innovative color coding and modular silicone holders:

- » The innovative shape reduces instrument contact allowing for more liquid flow and for more efficient cleaning.
- » 5 color coding tabs for practice organization included with every cassette.
- » Silicone rails can be assembled in any position, are compatible with all handle diameters.



1/2 DIN cassette for 10 instruments
PP40372-2



1/4 DIN cassette for 5 instruments
PP40371-2

Item #	Description
PP400371-2	Cassette for 5 instruments, 1/4 DIN, blue rails, incl. 5 color coding option (purple, red, yellow, green, orange, black), 185x80x37mm
PP400372-2	Cassette for 10 instruments, 1/2 DIN, blue rails, incl. 5 color coding option (purple, red, yellow, green, orange, black), 185x148x37mm
PP400373-3	Cassette for 20 instruments, DIN, blue rails, incl. 5 color coding option (purple, red, yellow, green, orange, black), 285x183x37mm

Cassettes will be delivered without instruments.

YOUNG INNOVATIONS



Young Innovations Europe GmbH
Kurfürsten-Anlage 1 | D-69115 Heidelberg | Germany
+353 (58) 45 966 | info@ydnt.eu

ORAL HYGIENE

PEDIATRIC

RESTORATIVE